



Cleaning and Sterilization for MPR Knee Positioning Triangles

Covers the entire range of Positioning Triangles from MPR

Recommended hand cleaning instructions

For the MPR-193 & MPR-293. Soak the product in an enzyme solution. Follow the manufacturer's directions for dilution and soaking time.

For the MPR-193P & MPR-293P. Please use your facilities recommended aluminum detergent, take care to follow the manufacturer's instructions.

Impacted material should be scrubbed out; however, steel wool or abrasive soap cleaners should not be used as they may scratch the surface finish of the device.

Rinse the product in warm tap water.

Soak or wipe the product down with a hospital approved and/or EPA approved germicide.

Rinse the product in warm tap water.

Wash the MPR-193 & MPR-293 with an instrument detergent up to a pH of 9.0. For the MPR-193P & MPR-293P. Please use your facilities recommended aluminum detergent, take care to follow the manufacturer's instructions.

Rinse the product in warm tap water.

Dry thoroughly and wrap.

Washer Cleaning Instructions

Place items through decontaminator /washer, per the manufacturer's instructions.

For the MPR-193P & MPR-293P. Please use your facilities recommended aluminum detergent, take care to follow the manufacturer's instructions.

Prepare items for sterilization.

Note: The device must be cleaned prior to sterilization.

Recommended Sterilizer Instructions

Sterilize components based upon your facilities protocol.

For the MPR-193P & MPR-293P. Please use your facilities recommended aluminum detergent, take care to follow the manufacturer's instructions.

Ensure that all parts are thoroughly cleaned.

Sterilizer Type: Pre-Vacuum Steam Sterilizer

Temperature: 132° C / 270° F
Minimum Cycle Time: 4 minutes
Minimum Drying Time: 30 minutes.

Can be pre-vacuum steam sterilized unwrapped or double wrapped in two disposable/reusable wraps. Run normal vacuum cycle for your institution.