


MPR Lateral Hip Positioning System set-up guide



Step 1:

The MPR Lateral Hip Positioning Peg Board System is placed on the surgical table (taking careful note to follow your facilities guidelines of either directly on the table, or on top of the table pad) and secured with the two stabilizing clamps provided — one at the head of the table, and the other at the foot of the positioner on the opposite side of the table. The clamps are secured and locked to the table rails using the MPR Universal Table Clamp, or any other blade style table clamp.

Step 2:

Place the gel pad on the positioner. A sheet (not shown) can then be placed over the padding to help prevent shifting as the patient is transferred to the table.

Step 3:

Pegs should be positioned as follows:

- 1) A short peg placed in the sacral area. A short peg is used here instead of a long peg in order to maximize the operative area.
- 2) A short peg placed in front of the pubis. Care should be taken to avoid neurovascular structure compression in both femoral triangles.
- 3) A long peg placed in the posterior thoracic region below the scapula.
- 4) A long peg placed directly below the pectoral muscle, keeping in mind the need for adequate ventilation of this area. *Note: More pegs, as well as the use of the peg extenders, may be required for larger patients.*

CAUTION: The board should not be used as a transport device.