

Vinyl-Covered Foam

Vinyl-covered foam products offer stain-resistant, anti-microbial, anti-bacterial, and waterproof properties, to aid in easier cleanability. It is important to note that this *does not* protect the patient from pathogens if proper procedures are not followed. Please refer to the following guidelines for cleaning:

1. Before cleaning, inspect the product thoroughly for any damage or tears in the vinyl. If any imperfections are found, remove the product from use until it has been repaired or replaced.
2. Thoroughly wipe the product down with a neutral detergent to remove any soil, then disinfect with a hospital- or EPA-approved germicide or quaternary solution according to your facility's protocols.
3. Allow to dry completely before use or storage. *Do not* use heat, or any harsh chemicals or abrasive materials; these may deteriorate the vinyl covering.

Closed-Cell Foam

Closed-cell foam products offer lightweight, durable, non-porous, antimicrobial, and antifungal properties without the need for a covering. It is important to note that this *does not* protect the patient from pathogens if proper procedures are not followed. Please refer to the following guidelines for cleaning:

1. Before cleaning, inspect the product thoroughly for any damage or tears. While damage may not affect product's usability, it is important to note. Consult your facility's protocols before use.
2. Thoroughly wipe the product down with a neutral detergent to remove any soil, then disinfect with a hospital- or EPA-approved germicide or quaternary solution according to your facility's protocols.
3. Allow to dry completely before use or storage. *Do not* use heat, or any harsh chemicals or abrasive materials; these may deteriorate the quality of the foam.

Spray-Coated Foam

Spray-coated foam products are stain-resistant and antimicrobial, aiding in easier cleanability. It is important to note that this *does not* protect the patient from pathogens if proper procedures are not followed. Please refer to the following guidelines for cleaning:

1. Before cleaning, inspect the product thoroughly for any damage or tears in the vinyl. If any imperfections are found, remove the product from use until it has been repaired or replaced.
2. Thoroughly wipe the product down with a neutral detergent to remove any soil, then disinfect with a hospital- or EPA-approved germicide or quaternary solution according to your facility's protocols.
3. Allow to dry completely before use or storage. *Do not* use heat, or any harsh chemicals or abrasive materials; these may deteriorate the quality of the coating.