



MPR-9020 Carter Hand Table System Instructions for Use

Set-Up and Use

It is crucial to conduct a thorough examination of the Carter Hand Table for any defects or damage. This is to ensure patient safety and the effectiveness of the table. Look for any visible signs of wear and tear, cracks, or loose components. Do not use if damage is observed.

Adhere to your facility's specific cleaning and sterilization protocols to decontaminate the table. This is a critical measure in maintaining a sterile environment.

On the underside of the table, there's a mounting lever. Make sure this lever is in an upward position, with its handle resting near the bottom surface of the table. This position ensures that the table is ready for mounting.

The Carter Hand Table is designed with locating slots designed for easy and secure attachment to the OR table. Align these slots with the rail on the OR table and slide the hand table into its designated position.

Once the table is positioned correctly, secure it by moving the mounting lever downwards, towards the floor. This action locks the Carter Hand Table in place, ensuring it remains stable during procedures.

Before use with the patient, it's important to test the setup. Apply a little pressure to the table to make sure it's secure and doesn't shift. This final check is crucial for ensuring the safety and comfort of the patient during the procedure.

Following these steps, you can ensure the safe and effective use of the Carter Hand Table.