



Cleaning and Sterilization - MPR Reusable Foam Products

Cleaning Guide for MPR reusable foam Products

Cleaning and Sterilization

MPR Reusable Foam products offer closed-cell foam construction.

Closed-cell foam allows easy cleaning and disinfection; the material is compatible with standard cleaning solutions.

Care should be taken to keep the reusable foam product clean by washing with a soapy solution and soft brush after each use.

After proper cleaning, store where it will dry easily.

Applies to the following MPR Foam Products

MPR-1503	MPR NS.0007
MPR-5501	MPR NS.0007C
10201	MPR NS.0009
10201C	MPR NS.0009C
10203	MPR NS.0011
10203C	MPR NS.0012
10204	MPR NS.0012-RWB
10205	MPR NS.0012-RWBXL
10233BC	MPR NS.0013
10225C	6823-14-000
10229C	6823-16-000
10233C	