



## MPR Golden Magnetic Retriever

### Description

A magnetic cannula attachment for removing metal fragments during procedures. Its slender design allows it to fit into any standard 4mm cannula and be easily manipulated within joints. The operative end houses a  $\frac{1}{4}$ " magnet, and the opposite end is open for the purpose of attaching a length of suction tubing. The gold finish makes it stand out among stainless steel instruments, ensuring quick identification in the event of a broken blade.

Made in USA. Designed by Dr. Lanny Johnson, M.D.

### Product Specifications

Overall Length	9"/22.86cm
Diameter	0.16"/4mm



### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

