



MPR Kool Grip Mallet



Description

An orthopedic mallet designed with reduced mass and material, minimizing hand fatigue during procedures and heat retention during sterilization. The interchangeable heads are available in either Delrin® or a special steel alloy, and can be purchased with the mallet handle, or separately (*see below for more information*).

MPR-6809	Kool Grip Mallet (handle only)
MPR-6809-M	Kool Grip Mallet w/ Metal Heads
MPR-6809-P	Kool Grip Mallet w/ Delrin® Heads

MPR-6809-M MPR-6809-P



REPLACEMENT HEADS

MPR-6809-MR	Set of Metal Heads (2 ct.)
MPR-6809-PR	Set of Delrin® Heads (2 ct.)

Product Specifications

Overall Length	9"/22.86cm
Handle Length	6.5"/16.51cm
Width (w/o Heads)	2.9"/7.41cm
Width (w/ Heads)	4"/10.31cm
Head Diameter	1.5"/3.81cm
Weight (w/o Heads)	1.05 lbs/476g
Weight (w/ Metal Heads)	1.85 lbs/839g
Weight (w/ Delrin® Heads)	1.15 lbs/522g

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.