



## MPR MIS LIT Offset Double Point Retractors



MPR-7804-103



MPR-7804-104

### Description

Retractors used during total hip replacement procedures using the Anterolateral Approach technique to retract the femur posteriorly. The design features a 90° curved blade, a double-pronged tip for a secure grip with square points to protect the patient, and two holes in the handle for attaching weights. A light pipe can also be attached via a groove on the underside and a hole in the blade.

All MPR Orthopedics retractors are made in the USA from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

### Product Specifications

#	DESCRIPTION	LENGTH	BLADE L.	BLADE W.	SHORT TIP*	LONG TIP*
MPR-7804-103	Left Twisted	14.97" (38cm)	5.51" (14cm)	1.25" (32mm)	0.71" (18mm)	1.08" (27.5mm)
MPR-7804-104	Right Twisted	14.97" (38cm)	5.51" (14cm)	1.25" (32mm)	0.71" (18mm)	1.08" (27.5mm)

\*measured from end of blade

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).