

MPR Bone-Holding Forceps, Narrow Jaw

MPR399.051



Description

Reduction forceps for holding and manipulating bones during surgery, ensuring precise alignment and stability. Their time-proven design features serrated jaws for a secure grip, and a speed-lock mechanism allowing them to easily be held in the desired position. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Comparable to Synthes Bone-Holding Reduction Forceps* (399.05 or 399.051).

Made in Germany.

Product Specifications

Overall Length	9.45" (24cm)
Tip Spread	1.69" (43mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Synthes is a registered trademark of that corporation.