

MPR “Captain Hook” Wire Twister

MPR 1292-97



Description

An instrument for quickly and easily twisting surgical wires. The freely rotating hook is made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and the Radel®* handle is ergonomically designed to provide a secure, comfortable grip.

Used with Luque wires, single wires with eyelets, and cerclage wires.

Made in USA.

Product Specifications

Overall Length	7.5" (19.05cm)
Handle Length	4.25" (10.8cm)
Handle Diameter	1.17" (2.97cm)
Hook Length	3.5" (8.89cm)
Hook Diameter	0.31" (8mm)
Tip Width	0.12" (3mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.