



MPR Jaws Graspers



10027



10027R



10028



10028R

Description

Graspers with cupped jaws for removing large loose bodies during procedures. Their low-profile design features elliptical barrels and tips for maximum visibility and precise control. MPR Orthopedics offers four different models, in “Grande” and “Venti” sizes, and with or without ratchet-lock mechanisms. The “Grande”-size models feature fine serrations on both jaws, and the “Venti”-size models feature five retrograde serrated teeth on the upper jaw. All models are designed to ensure a secure grip and prevent tissue from watermeloning. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

#	SIZE	LENGTH	RATCHET-LOCK	TIP LENGTH	TIP WIDTH	TIP SPREAD
10027	Grande	6.75" (17.15cm)	No	0.47" (12mm)	0.16" (4mm)	0.41" (10.5mm)
10027R	Grande	6.75" (17.15cm)	Yes	0.47" (12mm)	0.16" (4mm)	0.41" (10.5mm)
10028	Venti	6.75" (17.15cm)	No	0.59" (15mm)	0.22" (5.6mm)	0.48" (12.25mm)
10028R	Venti	6.75" (17.15cm)	Yes	0.59" (15mm)	0.22" (5.6mm)	0.48" (12.25mm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.