

MPR Straight Rasp

10705



Description

A rasp for restimulating the blood supply for complete meniscal repair in both open and arthroscopic procedures. The knurled grip allows for precise control, and its unique tooth design allows it to cut while pushing and pulling. Also used in conjunction with a suture passer in the Meniscus Mender II system.

Made in USA. Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	9" (22.86cm)
Diameter	0.12" (3mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.