



MPR Phenolic Handle Mallets



21168

21174

Description

Universal orthopedic mallets. The grooved handles are made from a durable, heat- and wear-resistant phenolic resin, and ergonomically designed to provide a secure, comfortable grip. The heads are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

MPR Orthopedics offers two models with different sizes and weights to suit your needs and/or preference. The handles can be safely steam-sterilized up to 290°F/143°C.

Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	WEIGHT	HEAD DIAMETER	HEAD WEIGHT
21168	MPR Phenolic Handle Mallet, 10"	10" (25.4cm)	27oz (765g)	1.38" (3.5cm)	16oz (455g)
21174	MPR Phenolic Handle Mallet, 11"	11" (28cm)	41oz (1160g)	1.77" (4.5cm)	32oz (900g)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.