

## MPR Beckman-Weitlaner Retractor

MPR-2557



### Description

A self-retaining retractor for spreading tissue to access a surgical site. The arms are hinged to allow use in anatomically challenging areas and the pronged tips ensure a secure grip on tissue.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

### Product Specifications

Overall Length	6.38" (27.5cm)
Arm Length (Shorter)	2.5" (6.35cm)
Arm Length (Longer)	2.63" (6.67cm)
Tip Length (Shorter)	0.5" (1.27cm)
Tip Length (Longer)	0.75" (1.91cm)
Prong Length	0.55" (14mm)
Prong Width	0.1" (2.5mm)
Tip Spread	3.88" (9.86cm)
Range of Rotation*	90°

\*both directions

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).