



MPR Offset Retractors



MPR-31-555514



MPR-31-555515

Description

Retractors with offset blades for use in orthopedic procedures. Other features include a 75° bent design, a rounded tip, and a hole at the bottom of the handle for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Product Specifications

#	DESCRIPTION	LENGTH	BLADE L.*	BLADE W.	TIP L.	TIP W.
MPR-31-555514	Offset Right	10.25" (25.4cm)	2.36" (6cm)	0.51" (13mm)	0.83" (21.13mm)	0.31" (8mm)
MPR-31-555515	Offset Left	10.25" (25.4cm)	2.36" (6cm)	0.51" (13mm)	0.83" (21.13mm)	0.31" (8mm)

*measured along outer edge

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.