

## MPR Patella Forceps

MPR-398-76



### Description

Forceps for holding and aligning the patella. The self-centering design features forked tips for a precise and secure grip, and a speed-lock mechanism for holding the instrument in the desired position. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Comparable to Synthes Patella Forceps\* (398.76).

### Product Specifications

Overall Length	7.13" (17.5cm)
Tip Spread	0.05"-2.01" (1.25mm-51mm)
Fork Width	0.33" (8.5mm)
Prong Width	0.03" (0.75mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

\*Synthes is a registered trademark of that corporation.