



## MPR Mini Hohmann Retractors



MPR-4101

MPR-4102

### Description

Miniature Hohmann-style retractors, specialized for delicate tissues. The tapered handle design and narrow 15°-bent blades maximize precision while minimizing visual obstruction, and the square points help to prevent tissue damage.

MPR Orthopedics offers two models with differing blade widths, to suit your needs or preferences. All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

### Product Specifications

#	DESCRIPTION	LENGTH	HANDLE	BLADE L.	BLADE W.	POINT
MPR-4101	6mm x 6.5"	6.3" (16cm)	4.5" (11.4cm)	1.14" (29mm)	0.24" (6mm)	0.12" (3mm)
MPR-4102	8mm x 6.5"	6.3" (16cm)	4.5" (11.4cm)	1.14" (29mm)	0.31" (8mm)	0.12" (3mm)

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).