

MPR Universal Guide Pin Hand Grip

Description

A grip for inserting guide wires into fractured bones. The highly versatile design features an ergonomic handle for a secure, comfortable grip, and can be easily disassembled for cleaning. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	3.75" (9.53cm)
Handle Width	3.5" (8.9cm)
Handle Diameter	0.63" (16mm)
Top Hole Diameter	0.2" (5mm)
Clamping Diameter	0.1"-0.14" (2.5mm-3.5mm)

MPR-469675



For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.