

## MPR Blunt Hohmann Retractor

MPR-6079-950



### Description

A specialized Hohmann retractor with a blunt, rounded tip and slender leaf-shaped blade. Its curved design provides a high degree of access and leverage, and the handle features three holes for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

### Product Specifications

|                |                  |
|----------------|------------------|
| Overall Length | 10.75" (27.31cm) |
| Handle Length  | 5.1" (12.95cm)   |
| Blade Length   | 2.9" (7.37cm)    |
| Blade Width*   | 0.89" (2.25cm)   |
| Tip Length     | 0.65" (1.65cm)   |
| Tip Width      | 0.38" (0.96cm)   |

\*measured at widest point

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).