

MPR Sharp Hohmann Retractor

MPR-6079-955



Description

A specialized Hohmann retractor with a sharp, pointed tip for fixing to bone. The narrow blade and tip are designed to minimize visual obstruction during procedures, and the handle features a pair of holes for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Product Specifications

Overall Length	9" (22.86cm)
Handle Length	6.5" (16.5cm)
Blade Length	1.38" (3.5cm)
Blade Width	0.67" (1.7cm)
Tip Length	1" (2.54cm)
Tip Width	0.31" (0.8cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.