

## MPR Kool Grip Mallet



### Description

An orthopedic mallet designed for minimal hand fatigue during procedures and heat retention during sterilization. The interchangeable heads are available in either Delrin® or a special steel alloy, and can be purchased with the mallet handle, or separately (*see below*).

Made in USA.

MPR-6809	Kool Grip Mallet (handle only)
MPR-6809-M	Kool Grip Mallet w/ Metal Heads
MPR-6809-P	Kool Grip Mallet w/ Delrin® Heads

### REPLACEMENT HEADS

MPR-6809-MR	Set of Metal Heads (2 ct.)
MPR-6809-PR	Set of Delrin® Heads (2 ct.)

MPR-6809-M MPR-6809-P



### Product Specifications

Overall Length	9" (22.86cm)
Handle Length	6.5" (16.51cm)
Width (w/o Heads)	2.9" (7.41cm)
Width (w/ Heads)	4" (10.31cm)
Head Diameter	1.5" (3.81cm)
Weight (w/o Heads)	1.05 lbs (476g)
Weight (w/ Metal Heads)	1.85 lbs (839g)
Weight (w/ Delrin® Heads)	1.15 lbs (522g)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).