

MPR Mikhail Collateral Ligament Retractor

MPR-86-9054



Description

A retractor used during knee-related procedures to protect the medial and collateral ligaments. Its slender, single-pronged tip is designed to maximize precision and access, with a rounded point, and the handle features a pair of holes for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

Product Specifications

| | |
|----------------|-----------------|
| Overall Length | 9.75" (24.77cm) |
| Handle Length | 6.75" (17.15cm) |
| Tip Length | 2.5" (6.35cm) |
| Tip Width* | 0.55" (1.4cm) |
| Point Width | 0.2" (0.5cm) |

*measured at widest point

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.