

MPR Tibial Retractor/PCL Retractor

MPR-9189



Description

A retractor used during total knee implant surgeries to expose the distal femur and proximal tibia to allow precise saw cuts, while protecting the posterior cruciate and medial/lateral ligaments. Features a dual-pronged circular tip, and a pair of holes in the handle for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

Product Specifications

Overall Length	10" (25.4cm)
Handle Length	7" (17.78cm)
Tip Length	1.25" (3.18cm)
Tip Width	1.25" (3.18cm)
Tip Spread	0.44" (1.12cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.