



# MPR Beach Chair Replacement Straps

MPR-F-STRAP



## Description

Replacement straps for Smith & Nephew's T-MAX Beach Chair Positioner\*, with buckles on either end and velcro on one side. The straps attach to the surgical headset of the chair via the velcro, and clip to the buckles on the face mask.

Made in USA. Sold in pairs.

## Product Specifications

Overall Length	11.5" (29.21cm)
Strap Width	1" (2.54cm)

## For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

\*Smith & Nephew and T-MAX are registered trademarks of Smith & Nephew.