




MPR Capsule Retractors



MPR1054-01

MPR1054-02

Description

Retractors used during hip orthopedic surgeries to retract the acetabular capsule and expose the acetabulum. Features include rounded tips for easy tissue retraction while protecting the patient, and a hole in the handle for attaching weights. MPR Orthopedics offers two different capsule retractors, for each hip.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

Product Specifications

#	DESC.	LENGTH	BLADE LENGTH	BLADE WIDTH	WIDE TIP WIDTH	THIN TIP WIDTH
MPR1054-01	Left	12" (30.5cm)	1.38" (3.5cm)	0.79" (2cm)	0.35" (0.9cm)	0.18" (4.5mm)
MPR1054-02	Right	12" (30.5cm)	1.38" (3.5cm)	0.79" (2cm)	0.35" (0.9cm)	0.18" (4.5mm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.