



MPR Single-Use Anti-Fatigue Pads



MPR106-171317-4



MPR108-171830-4

Description

Single-use anti-fatigue mats to provide support and reduce discomfort/fatigue from standing for long periods of time. Their durable polyurethane surface resists fluids, dirt, and staining, and is textured to provide friction, and the beveled edges minimize risk of tripping. The one-time use design reduces risk of contamination and greatly improves OR turn-over time.

MPR Orthopedics offers two different styles of single-use anti-fatigue mats, a floor mat and a stool mat. The latter is intended to be used in conjunction with a step stool such as the MPR Stackable Step Stool (MPR140; pictured, not included). These mats are intended for single use *only*; discard after use.

Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	WIDTH	THICKNESS	COUNT
MPR106-171317-4	Single-Use Floor Mat	18" (45.72cm)	30" (76.2cm)	0.63" (1.58cm)	5
MPR108-171830-4	Single-Use Stool Mat	13" (33.02cm)	17" (43.18cm)	0.63" (1.58cm)	10

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.