

MPR Medial Calcar Retractor

MPR110031569



Description

A retractor used during anterior hip procedures such as anterior hip arthroplasty and minimally invasive procedures to expose the medial calcar. The balanced design ensures a comfortable grip and precise control, and features a slender blade with a dual-pronged tip for stable retraction with minimal obstruction, and a pair of holes in the handle for attaching weights.

Designed to facilitate the Anterior Supine Inter-muscular (ASI) approach.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Product Specifications

Overall Length	15" (38.2cm)
Handle Length	9" (23.5cm)
Blade Length	3.34" (8.5cm)
Tip Width	1.18" (2.99cm)
Bend	25°

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.