



# MPR Stainless Steel Anesthesia Tube Tree

MPR110



## Description

A tube tree for routing tubes, lines, cords, and cables between a patient and medical equipment. Designed to be tucked under a patient or surgical bed mattress, it eliminates the need for taping or pinning lines while keeping the space neat and lines organized. Features two sides with six and four spaces for tubes respectively, depending on the needs of the operation and/or space.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

## Product Specifications

Length (6-space side)	13" (33.02cm)
Length (4-space side)	8.5" (21.59cm)
Width	5" (12.7cm)

## For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).