



MPR Stainless Steel Mallets, Classic Design



MPR155-01



MPR155-02



MPR155-05

Description

Universal orthopedic mallets. Their time-proven design has been sold for decades by leading orthopedic manufacturers, and is optimized for a balanced feel, with a knurled handle to provide a secure, comfortable grip.

MPR Orthopedics offers three different models with a range of sizes and weights. All are made from solid stainless steel, ensuring durability and longevity, and allowing for easy cleaning/sterilization.

Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	HANDLE	HEAD Ø	WEIGHT
MPR155-01	Stainless Steel Mallet, 1 lb.	7.63" (19.37cm)	4.5" (11.4cm)	1.25" (3.18cm)	1 lb (454g)
MPR155-02	Stainless Steel Mallet, 2 lbs.	7.75" (19.69cm)	4.5" (11.4cm)	1.5" (3.81cm)	2 lb (907g)
MPR155-05	Stainless Steel Mallet, 5 lbs.	10.25" (26.35cm)	6.5" (16.5cm)	2.5" (6.35cm)	5 lb (2,268g)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.