

# MPR Paulson Blount Knee Retractor

MPR280800000



## Description

A retractor for use during knee procedures such as meniscectomies and arthrotomies. The ring-shaped handle minimizes heat retention after sterilization, and the thin, shallowly curved blade maximizes precision and dexterity for use on more delicate tissues.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

## Product Specifications

Overall Length	7.28" (18.5cm)
Handle Length	5.24" (13.3cm)
Handle Width	1.25" (3.18cm)
Handle Diameter	0.31" (0.8cm)
Blade Length	2.77" (7.03cm)
Blade Width	0.51" (1.3cm)

## For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).