



MPR Closed Case Carts



Description

Carts for safe, sterile transportation of surgical instruments between Sterile Processing and the OR. Feature integrated handles, adjustable internal shelves, directional- and total-lock wheel casters, and door catches to make cleaning in cart washers easier.

MPR Orthopedics offers single- and double-door models, with either one or two internal shelves. Designed and made from materials proven to withstand the conditions of cart washers and surgical settings.

Made in USA.

Product Specifications

SINGLE-DOOR MODELS

#	EXT. HEIGHT	EXT. WIDTH	EXT. DEPTH	INT. HEIGHT	INT. WIDTH	INT. DEPTH
MPR3040	39" (99.06cm)	37" (93.98cm)	27" (68.58cm)	27" (68.58cm)	27.5" (69.85cm)	24" (60.96cm)
MPR3045*	59" (149.86cm)	37" (93.98cm)	27" (68.58cm)	47" (119.38cm)	27.5" (69.85cm)	24" (60.96cm)

DOUBLE-DOOR MODELS

#	EXT. HEIGHT	EXT. WIDTH	EXT. DEPTH	INT. HEIGHT	INT. WIDTH	INT. DEPTH
MPR3050	39" (99.06cm)	48" (121.92cm)	27" (68.58cm)	27" (68.58cm)	38.5" (97.79cm)	24" (60.96cm)
MPR3055*	59" (149.86cm)	48" (121.92cm)	27" (68.58cm)	47" (119.38cm)	38.5" (97.79cm)	24" (60.96cm)

*two shelves

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.