

MPR Dual Z-Retractors



MPR3313



MPR3313-32

Description

Retractors for exposing the femur and tibia during knee-related procedures, such as joint replacement. The distinctive low-profile “Z”-shaped design offers a high degree of leverage while keeping the handle clear of the operative field and providing maximal access to the articulating surfaces. MPR Orthopedics offers two models, with different shaft lengths.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA

Product Specifications

#	DESCRIPTION	LENGTH	HANDLE L.	HANDLE W.	BLADE L.	BLADE W.
MPR3313	Dual Z Retractor, Long	7" (17.78cm)	3" (7.62cm)	0.79" (2cm)	2.5" (6.35cm)	0.52" (1.33cm)
MPR3313-32	Dual Z Retractor, Short	3.5" (8.9cm)	3" (7.62cm)	0.79" (2cm)	2.5" (6.35cm)	0.52" (1.33cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.