

MPR Minor Procedure Arm and Hand Surgery Table

MPR333



Description

A surgical table extension for arm and hand procedures. Its sturdy phenolic construction ensures durability and longevity, and its low-profile legless design attaches directly to a US surgical table rail via a pair of built-in clamps, making it perfect for C-Arm imaging.

Comes with a 2"-thick reusable pad (MPR334, *second picture*).

Made in USA.

Product Specifications

Length	24" (60.96cm)
Width	15" (38.1cm)

Accessories

#	DESCRIPTION
MPR334	2" Replacement Pad
MPR334-4	4" Replacement Pad



For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.