

## MPR Periosteal Elevator

MPR399.40



### Description

An instrument for separating periosteal tissue from bone during fracture repair and other such procedures. The low-profile design features a straight-edged curved blade, ensuring precise and clean tissue separation, and a faceted Radel®\* handle to prevent rolling when placed down.

Made in USA.

### Product Specifications

|                 |                 |
|-----------------|-----------------|
| Overall Length  | 7.31" (18.58cm) |
| Handle Length   | 3.93" (10cm)    |
| Handle Diameter | 0.63" (16mm)    |
| Shaft Length    | 2.5" (6.35cm)   |
| Shaft Diameter  | 0.2" (5mm)      |
| Blade Length    | 1" (2.54cm)     |
| Blade Width     | 0.25" (6mm)     |

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

\*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.