

MPR Femoral Extractor Slap Hammer



Description

A device used in conjunction with an extractor for implant stem removal. The sliding weight is ergonomically designed to ensure a secure, comfortable grip, and the spring loaded sleeve around the head allows for easy attachment to the selected extractor, and holds it firmly in place. Made from solid stainless steel for longevity, durability, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	19" (48.26cm)
Handle Length	3.5" (8.89cm)
Weight Length	6.63" (16.84cm)
Overall Weight	3.65 lbs. (1.66kg)

Product Compatibility

#	DESCRIPTION
MPR944	MPR Universal Femoral Stem Extractor
MPR945	MPR One-Piece Stem Adaptor
MPR946	MPR Closed Loop Extractor
MPR947	MPR Hook Stem Extractor

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.