

MPR Ortho Rongeur with Pistol Grip

05523



Description

A rongeur for cutting and shaping of bone during revision procedures. The pistol-grip design ensures stable handling and precision and is engineered for smooth gripping action, while the broad-tipped jaws and double-jointed clamping mechanism maximize cutting strength. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	8.75" (22.23cm)
Grip Length	3.25" (8.26cm)
Handle Width	0.38" (0.95cm)
Tip Length	0.63" (1.6cm)
Tip Width	0.31" (8mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.