



## MPR Regular Probe

10014



### Description

A standard orthopedic probe. The knurled handle ensures a secure, comfortable grip, and the larger tip compared to other probes makes it the ideal probe for all purposes. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

### Product Specifications

Overall Length	9" (22.86cm)
Shaft Diameter	0.13" (3.2mm)
Tip Width	0.12" (3mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).