




## MPR Single-Use Carpal Tunnel Knives



10030

10030-15

### Description

Single-use knives for cutting into the wrist to access the carpal tunnel. The razor-sharp blades are set at a 65° angle to eliminate the need for any sawing motion that may lead to loss of control and double-shielded to prevent any unintended damage to the nerves and tendons. Two different models are available, with either a straight or 15°-bent design. Used in conjunction with the MPR Carpal Tunnel Retractor (10031; *sold separately*).

Intended for single use *only*; safely dispose after use.

Made in USA. Sold in packs of 10.

### Product Specifications

#	DESCRIPTION	OVERALL LENGTH	HANDLE LENGTH	BLADE LENGTH	BLADE WIDTH
10300	Straight	7.5" (19.05cm)	4.8" (12.2cm)	0.2" (5mm)	0.33" (8.5mm)
10300-15	Bent	7.5" (19.05cm)	4.8" (12.2cm)	0.2" (5mm)	0.33" (8.5mm)

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).