

MPR Lightweight Mallet

21133



Description

An orthopedic mallet designed to be as light-weight as possible without sacrificing dependability or usability. The handle is constructed from anodized aluminum, and the head is solid German stainless steel. Comes with a pair of replaceable nylon caps.

Made in USA.

Product Specifications

Overall Length	7.5" (19.1cm)
Overall Weight	7oz (200g)
Handle Length	3.54" (9cm)
Handle Width	0.79" (2cm)
Head Width (base)	1.57" (4cm)
Head Width (w/ caps)	2.36" (6cm)
Head Diameter	0.98" (2.5cm)
Head Weight	4oz (120g)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.