

MPR Dead Blow Mallet

21162



Description

A repercussion-free orthopedic mallet for revision procedures. Features a grooved phenolic resin handle for a secure, comfortable grip, and a head with one nylon cap* and one stainless steel face. Made from stainless steel for durability, longevity, and easy cleaning/sterilization. The handle can be safely steam-sterilized up to 250°F/121°C.

Made in USA.

Product Specifications

Overall Length	11.5" (29.21cm)
Overall Weight	21oz (600g)
Handle Length	5.25" (13.34cm)
Handle Diameter	1.13" (2.86cm)
Head Width	3.5" (8.89cm)
Head Diameter	1.26" (3.2cm)
Head Weight	10oz (285g)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Nylon cap is not autoclavable; remove before sterilizing.