

MPR Four Prong Rake Retractor

MPR-2810-20



Description

A retractor used for pulling and holding open the edges of wounds and incisions. The deeply curved rake-shaped tip ensures a secure grip on skin and subcutaneous tissue, while the long, slender shaft minimizes clutter around the operative space. The exclamation-point shaped handle is designed to be as ergonomic and light-weight as possible, and for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

Product Specifications

Overall Length	8.66" (22cm)
Handle Length	4.44" (11.28cm)
Tip Length	1" (2.54cm)
Tip Width	0.75" (1.91cm)
Prong Width	0.16" (4mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.