

## Modified Stille-Luer Rat Teeth Bone Rongeur

MPR-3314



### Description

A modified version of a Stille-Luer style bone rongeur, for gouging out and removing smaller bone fragments during procedures. The “rat-tooth” serrated tips are designed to dig into bone tissue with ease, and the double-leaf-spring mechanism ensures smooth and precise gripping action. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	8.5" (21.59cm)
Handle Width	4.25" (10.8cm)
Tip Width	0.31" (8mm)
Tip Thickness	0.22" (5.5mm)
Tip Spread	0.77" (2cm)
Tooth Depth	0.04" (1mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).