

## MPR Reduction Forceps, Angled Points

MPR-378



### Description

Self-retaining forceps used in fracture repair procedures to reposition and realign fractured bones. The sharp inward-turned points dig into bone tissue, ensuring a firm, secure grip, and the speed-lock mechanism allows it to be easily retained in and released from the desired position.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

### Product Specifications

Overall Length	8" (20.32cm)
Arm Length	2.75" (6.99cm)
Arm Diameter*	0.2" (0.51cm)
Tip Length	0.39" (1cm)
Max. Point Spread	3" (7.62cm)

\*measured at midpoint

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).