

MPR Femoral Lifter/Distractor

MPR-4844



Description

Also known as the Mikhail U-Shaped Femoral Distractor; an instrument used in hip procedures to distract the femur from the tibia. The knurled handle ensures a secure, comfortable grip, and the slender shaft easily slides between the femur and tibia, and is designed to provide optimal leverage when distracting the bone.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9.75" (24.77cm)
Overall Width	5.78" (14.68cm)
Handle Length	5.06" (12.85cm)
Handle Diameter	0.95" (2.41cm)
Shaft Length (short end)*	6.38" (16.21cm)
Shaft Diameter	0.31" (8mm)

*measured from butt of handle

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.