




MPR Femoral Head Gauges



MPR-6791-7-101



MPR-6791-7-102



MPR-6791-7-103

Description

Gauges for measuring the femoral head during hip arthroplasty procedures. All are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

#	LENGTH	WIDTH	THICKNESS	MIN. HOLE Ø	MAX. HOLE Ø
MPR-6791-7-101	5.75" (14.61cm)	7.75" (19.69cm)	0.12" (3mm)	1.5" (38mm)	1.81" (46mm)
MPR-6791-7-102	5.75" (14.61cm)	7.75" (19.69cm)	0.12" (3mm)	1.85" (47mm)	2.05" (52mm)
MPR-6791-7-103	5.75" (14.61cm)	7.75" (19.69cm)	0.12" (3mm)	2.09" (53mm)	2.48" (63mm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.