

## MPR Straight Ball Spike

MPR-7117



### Description

A straight ball spike for reducing and holding fractured bones in place. The sharp point is designed to dig into bone tissue, eliminating the risk of slippage, and the faceted Radel®\* handle prevents rolling when placed down. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	11.81" (30cm)
Handle Length	5.51" (14cm)
Handle Diameter	1" (2.54cm)
Shaft Length	5.1" (12.95cm)
Shaft Diameter	0.31" (0.8cm)
Ball Diameter	0.31" (0.79cm)
Tip Length	0.1" (0.25cm)
Tip Diameter†	0.1" (0.23cm)

†measured at base

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

\*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.