

## MPR Ortho Hand Starter Rasp

MPR-7803-65



### Description

A rasp for use in implant alignment to shape the medullary canal, smooth bone surfaces, and initiate canal entry. The handle features two finger-holes for precise control and the curved rasping tip features a rounded end for smooth entry into the medullary canal. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	12.64" (32.1cm)
Handle Length	3.3" (8.37cm)
Handle Width	1" (2.54cm)
Shaft Length	3.71" (9.42cm)
Shaft Width	0.2" (5mm)
Tip Length	3.98" (10.1cm)
Tip Width	0.2" (5mm)
Thickness	0.28" (7mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).