




MPR C-Thru Armboards



MPR-9026

MPR-9039

Description

Transparent armboards that slip underneath the pad of any stretcher or surgical table. The radioluscent acrylic construction and simple obstruction-minimizing design make them ideal for C-Arm imaging and easy to clean and disinfect.

MPR Orthopedics offers C-Thru Armboards in two different sizes, which can be used with or without the corresponding reusable pads (*sold separately, see ACCESSORIES*). Made in USA.

Product Specifications

#	DESCRIPTION	SHORT END LENGTH	LONG END LENGTH	WIDTH
MPR-9026	C-Thru Armboard	20" (50.8cm)	25" (63.5cm)	11" (27.94cm)
MPR-9039	C-Thru Hand Table	20" (50.8cm)	25" (63.5cm)	15.5" (39.37cm)

Accessories

#	DESCRIPTION	LENGTH	WIDTH	THICKNESS
MPR-9034	C-Thru Armboard Pad	20.5" (52.07cm)	11" (27.94cm)	2" (5.08cm)
MPR-9040	C-Thru Hand Table Pad	20" (50.8cm)	15.5" (39.37cm)	2" (5.08cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.