



MPR K-Wire Cases



MPR-2071-10



MPR-424-80



MPR-424-85



MPR-424-90

Description

Folding cases/racks for holding Kirschner Wires, commonly known as “K-Wires”, and Steinmann Pins. Improved from an older design, these cases offer up to more than three times the carrying capacity, while being optimized for compactness and convenience, and are labeled with both standard and metric units.

MPR Orthopedics offers three different case sizes, as well as a miniature rack with a convenient tilt-stand design. All are compatible with all manufacturers’ K-Wires and Steinmann Pins, and are made from solid stainless steel to ensure durability and longevity.

Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	WIDTH	HEIGHT*
MPR-2071-10	Mini K-Wire Rack	---	4.33" (11cm)	4.75" (12.07cm) [†]
MPR-424-80	Mini K-Wire Case	4.5" (11.43cm)	3.63" (9.22cm)	2.5" (6.35cm)
MPR-424-85	K-Wire Case	6.5" (16.5cm)	3.63" (9.22cm)	4.13" (10.49cm)
MPR-424-90	Large Pin/Wire Case	9.63" (24.46cm)	6.19" (15.72cm)	5.5" (13.97cm)

*open/unfolded
[†]including handle

See reverse for information on carrying capacities.

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

Carrying Capacities

MPR Mini K-Wire Rack

MPR-2071-10

The MPR Mini K-Wire Rack is designed to hold 4" (10.16cm) wires in the following quantities:

DIAMETER	#
0.028" (0.7mm)	42
0.035" (0.9mm)	39
0.045" (1.1mm)	48
0.054" (1.4mm)	24
0.062" (1.6mm)	22
TOTAL:	175

MPR Mini K-Wire Case

MPR-424-80

The MPR Mini K-Wire Case is designed to hold 3" (7.62cm) and 4" (10.16cm) wires in the following quantities:

DIAMETER	#
0.028" (0.7mm)	19
0.035" (0.9mm)	21
0.045" (1.1mm)	25
0.054" (1.4mm)	15
0.062" (1.6mm)	25
TOTAL:	105

MPR K-Wire Case

MPR-424-85

The MPR K-Wire Case is designed to hold 5" (12.7cm) and 6" (15.24cm) wires/pins in the following quantities:

DIAMETER	#
0.028" (0.7mm)	30
0.035" (0.9mm)	28
0.045" (1.1mm)	36
0.054" (1.4mm)	30
0.062" (1.6mm)	31
0.078" (2.0mm)	29
TOTAL:	184

MPR Large K-Wire Case

MPR-424-90

The MPR Large K-Wire Case is designed to hold 9" (22.86cm) wires/pins in the following quantities:

DIAMETER	#
0.035" (0.9mm)	8
0.045" (1.1mm)	9
0.062" (1.6mm)	8
0.078" (2.0mm)	8
0.094" (2.4mm)	9
0.110" (2.8mm)	9
0.125" (3.2mm)	8
0.142" (3.6mm)	8
0.157" (4.0mm)	8
0.189" (4.8mm)	8
TOTAL:	83