




MPR Shoulder Traction Hooks



MPR-STH



MPR-STH-Z

Description

Hooks for attaching traction weights during shoulder-related procedures, designed to work with most shoulder distractors.

MPR Orthopedics offers large and compact models, depending on your preference. Both are made of stainless steel for durability and longevity, and designed to hold multiple traction weights. Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	WIDTH*	DIAMETER
MPR-STH	Shoulder Traction Hook	6.89" (17.5cm)	5.12" (13cm)	0.25" (6.35mm)
MPR-STH-Z	Shoulder Traction Hook, Compact	5.12" (13cm)	3.54" (9cm)	0.25" (6.35mm)

*measured at widest point

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.